

***“Train Smarter.  
Ride Stronger..”***

**BK-24**  
Sport Range Ergometer



The smart ergometer designed  
for wireless exercise testing



**COSMED**  
The Metabolic Company

“COSMED BK-24 was developed together with well-known cyclists coaches and scientist to get a product with a perfect cycling feel that needs to be experienced”

- ▶ Pendulum system with electronic control for guaranteed accuracy
- ▶ Compatible with all ANT+ and BLE devices with FE-C and FTMS protocols
- ▶ Integrated ANT+ heart rate signal receiver
- ▶ Step-less seat and bar adjustment for the best bike fit
- ▶ Electronic shifters and sprint buttons, for end-user configuration
- ▶ Calibrated Watts in display
- ▶ Load controlled with Watts (speed independent) or force (in kp, speed dependent)



BK-24 is a new performance bike fully integrated with COSMED ecosystem, for complete metabolic assessment for sports applications.

Especially developed to reproduce a realistic cycling feel, this bike allows you to ride in the lab as you would in the field.

Smart connectivity, granted by both BLE and ANT+ protocols (FTMS and FE-C profiles), creates new opportunities to connect and communicate with third-party devices and applications.

These characteristics, combined with the innovative display, electronic shifters and sprint buttons, make BK-24 the perfect solution for exercise testing, controlled training or virtual racing.

## Usability

The robust ergometer together with the characteristic flywheel creates a unique and smooth cycle feeling.

The electric servo enables third parties to take control of the resistance of the bike. Both the shifters and sprint buttons can be customized to best suit the rider's requirements.

## Accuracy

Accurate measurement is provided by the electronic control of our unique pendulum system, that allows simple calibration of the bike and reproducibility of measurements that is erroneously taken for granted with many high-power ergometer applications.

## Connection

The display communicates with both Bluetooth and ANT+, via the standardized communication protocols Fitness Machine Service (FTMS) and Fitness Equipment Controls (FE-C), respectively, so that third-party applications on Android (Samsung, Sony, Huawei and more), iOS (Apple), Windows and Mac are all supported.

## Customization

Use your own saddle, handlebars, pedals, or choose among different options; including a handlebar with gear shifters and sprint button, longer adjustable seat post, adjustable handlebar and mountain bike handlebars.



Height-adjustable seat post with rail notch indicators allows unlimited height settings





Unique pendulum system allows simple calibration of the bike by means of the kp scale and a calibration weight



Standard combi SPD pedals with toe clip cage design



Drop bar handlebars designed for a true racing position

Technical Specifications	
Frame	
Handlebar	Road racing handlebar
Stem	Adjustable horizontal 120 mm, vertical 500-910 mm
Seatpost	Adjustable horizontal 120 mm, vertical 500-910 mm
Saddle	Black, Racing
BB	68/107 mm
Crank	Steel, 172.5 mm
Pedals	9/16", combi Cage/ SPD
Resistance	
System	Pendulum system
Workload range	15-1400 W (@200 rpm)
Flywheel	20 kg
Display	
Parameters	RPM, HR (ANT+), time, speed, resistance, distance, Power/Watt
Connectivity	
Wired connection	USB
Wireless connection	BLE and ANT+
Weight and Measurement	
Max user weight	150 Kg
Length	1405-1630 mm
Width	640 mm
Height	1265 mm at display, 1240 mm at saddle
Weight	77 kg
Power requirements	
Voltage	100-240V (12V, 5,0 A)



Innovative graphic display for quick and precise adjustment of the resistance and wireless connectivity



#### ITALY - Headquarters

COSMED Srl  
Rome  
+39 06 931-5492  
[info@cosmed.com](mailto:info@cosmed.com)

#### ITALY

COSMED Srl  
Milan  
+39 02 99765-920  
[milano@cosmed.com](mailto:milano@cosmed.com)

#### GERMANY

COSMED Deutschland GmbH  
Schweinfurt  
+49 (0)9735 81390 00  
[DE@cosmed.com](mailto:DE@cosmed.com)

#### FRANCE

COSMED France SASU  
Brignais  
+33 (0)4 478628053  
[FR@cosmed.com](mailto:FR@cosmed.com)

#### THE NETHERLANDS

COSMED Benelux BV  
Nieuwegein  
+31 (0) 88 10 50 500  
[BNL@cosmed.com](mailto:BNL@cosmed.com)

#### DENMARK

COSMED Nordic ApS  
Odense  
+45 6595 9100  
[DK@cosmed.com](mailto:DK@cosmed.com)

#### SWITZERLAND

COSMED Switzerland GmbH  
Fehraltorf  
+41 (0)43 50 869 83  
[CH@cosmed.com](mailto:CH@cosmed.com)

#### USA

COSMED USA, Inc.  
Concord, Chicago  
+1 800 4263763 Toll Free  
[USA@cosmed.com](mailto:USA@cosmed.com)

#### AUSTRALIA

COSMED Asia-Pacific Pty Ltd  
Artarmon  
+61 449 971 170  
[ANZ@cosmed.com](mailto:ANZ@cosmed.com)

#### HONG KONG

COSMED HK Ltd  
Kowloon  
+852 3708 3126  
[HK@cosmed.com](mailto:HK@cosmed.com)

#### CHINA

COSMED Guangzhou Medical  
Technology Co. Ltd.  
Guangzhou  
[CHINA@cosmed.com](mailto:CHINA@cosmed.com)

© COSMED

E & OE. Subject to alterations without prior notice.

Products may not be available in your region as depending on countries and certifications.

All trademarks, registered trademarks and logos are the property of their respective owners.




#### COSMED Srl

Via dei Piani di Monte Savello 37  
Albano Laziale - Rome 00041  
Italy  
+39 (06) 931-5492 Phone  
+39 (06) 931-4580 Fax

[cosmed.com](http://cosmed.com)

#### Distributed by



To know more

